


The new state-of-the-art MicroFlo™ capnograph specifically developed for veterinarians!

MicroFlo™ Capnograph

CAPNO-5



MICROFLO™ SIDE-STREAM CAPNOGRAPHY



Developed for point of care capnography monitoring of intubated or non-intubated patients in surgery, dentistry, procedural sedation, pre-op preparation, anesthesia recovery, and intensive / critical care, the Capno-5 provides:

- Real-time CO₂ waveforms (capnogram)
- End-tidal CO₂ (ETCO₂)
- Fractional inspired CO₂ (FiCO₂)
- Respiratory rate (RR)
- Audible and visual apnea alarms
- Delivered with table mount



With the very low sample flow rate of 50 ml per minute, the Capno-5 is suitable for use in a wide range of veterinary patients, including cats and dogs of all sizes!

And unlike all mainstream capnographs, which are not designed for use in dentistry and non-intubated patients, the Capno-5 can be used in all procedures for intubated or non-intubated patients.

Technical features:

The Capno5 includes:

- Moisture trap
- Airway Adapter
- AC power adapter
- Table mount

Data ranges

- Sample flow rate: 50 ml / min.
- Respiration range: 2 -150 bpm
- Dimensions: 155 x 68 x 35 mm
- Weight: 340 gr
- Display: 3.5'' (8.9 cm)
- Battery autonomy: up to 10 hours



- ✓ EASY-TO-USE CAPNOGRAPH
- ✓ LONG-LIFE BATTERY (UP TO 10 HOURS OF USE)
- ✓ LARGE DISPLAY 3.5'' (8.9 CM)



 **Mano Médical**
Medical equipment for veterinarians

22 rue Charles de Gaulle
22100 Lanvallay - France
Tel. +33 (0)2 96 85 86 76

export@manomedical.com

www.manomedical.com