

# VET Ultrasound System

**Pocket Scan** / by  
Mano Médical

FROM POCKET  
TO HAND  
ULTRASOUND  
SYSTEM



Mano  
Médical  
by génia

# Wireless color doppler ultrasound scanner.



Functional and ergonomic use in the clinic and also ideal companion for visits outside and emergencies

## • Dual-probe technology

- ▶ Micro-convex probe
- ▶ Linear probe

## • Direct WIFI connection

- ▶ To your smartphone and your tablet
- ▶ Internal 5G WIFI network

## • Go anywhere diagnostic tool

- ▶ Up to 5 hours working time
- ▶ Less than 2 hours inductive recharge time

## • Handheld device

- ▶ Pocket size and light weight (only 250 gr)
- ▶ Packed in a carrying case

## • Imaging power

- ▶ Color doppler
- ▶ Pulsed doppler



Specifications:

▶ Frequencies	Micro-convex: from 3.5 to 5 MHz Linear: from 7.5 to 10 MHz
▶ Elements	128
▶ Radius and scan angle	Micro-convex: R20 mm / 88° Linear 40 mm
▶ Scan depth	90/130/160/200 mm & 20/40/60/100 mm
▶ Display modes	B, B/W, color, PW, PDI
▶ Measures	Distance, area, circumference...
▶ Working time	Up to 5 hours
▶ Reload mode	Inductive (wireless) in less than 2 hours
▶ Dimensions / weight	156 x 80 x 20 mm / 250 g
▶ Connectivity	Smartphone and tablet (InternalWi-Fi)
▶ WiFi	5 G
▶ Operating system	Apple iOS, Android, Windows



22 rue Charles de Gaulle  
22100 Lanvallay - France  
Tel. +33 (0)2 96 85 86 76

[export@manomedical.com](mailto:export@manomedical.com)

[www.manomedical.com](http://www.manomedical.com)