

# apAlert®

Portable apnoea detector

## Easy anesthesia respiration monitoring



The **apAlert® respiration monitor** makes care of veterinary patients much easier by signalling each breath, measuring the time between breaths and by alarming if breathing stops, slows or becomes shallower.

*Respiratory depression is a prime indicator of anesthetic depth which explains the **apAlert®** effectiveness.*

In essence, the **apAlert®** is measuring oxygen supply to the patient, effectively another way of doing what a pulse oximeter does without the fuss, bother and cost.

- Audible and visual respiratory beep
- Digital display of time between breaths
- No-breath alarm adjustment 10-60 secs
- Hands-free reset (optical)
- Sensitivity covers kittens to horses
- Adjustable volume
- Rechargeable battery (up to 70 hours per charge)
- Low battery alarm



**apAlert®** is supplied with a sterilizable sensor, a 9 volt transformer and a mounting bracket.

**They breath easy...you breath easy!**

 **Mano  
Médical**  
Medical equipment for veterinarians

22 rue Charles de Gaulle  
22100 Lanvallay - France  
Tél. +33 (0)2 96 85 86 76

[export@manomedical.com](mailto:export@manomedical.com)

[www.manomedical.com](http://www.manomedical.com)