

# DENTISTRY

## Veterinary Specific Dentistry Devices

### ▶ Strong Micromotor including polishing function

Type E micromotor - light, compact and handy for interventions on small animals or devitalized teeth



- Regular speed from 0 to 35 000 turns/min.
- 1:1 contra angle
- Setting drive
- Compact and light type E motor
- Including handpiece troy and pedal



### ▶ Vet 200 Plus Ultrasonic scaler

Piezo-electric scaler - portable and powerful, perfect for veterinary use!



- Titanium base alloy handpiece
- 5 piezo-ceramic tips
- Automatic frequency tracking digitally controlled
- Easy installation
- Compact and light design
- CE certified

Including tips:



 Mano Médical

Medical equipment for veterinarians

22 rue Charles de Gaulle  
22100 Lanvallay - France

Tél. +33 (0)2 96 85 86 76

[export@manomedical.com](mailto:export@manomedical.com)

[www.manomedical.com](http://www.manomedical.com)