

ULTRASONIC SCALER D7.

Including water tank and Led light

Detachable and autoclavable handpiece with light for better visibility

- ▶ Scaling, perio, endo
- ▶ Clinical solutions are applicable for automatic water supply mode, improving the performance of perio and endo function dramatically
- ▶ Less heat from tip during operation
- ▶ Improves comfort degree of patients

Tips:



GD1: Removes calculus plaque of supragingival, neck of teeth and hard calculus on teeth point



GD4: Removes all supragingival calculus and plaque



GD5: Removes all supragingival, neck calculus and plaque



GD6: Removes all supragingival heavy calculus



PD1: Removes subgingival calculus



PD3: Removes subgingival and deep perio calculus



PD4: Removes subgingival calculus and interdental calculus



Supplied with tip holder



Technical Data:

| | |
|-------------------|--|
| Power input: | 110V~ 50Hz/60Hz 220 - 240V~ 50Hz/60Hz |
| Main unit input: | 24V~ 50Hz/60Hz 1,3A |
| Output power: | 3W~20W |
| Frequency: | 28kHz \pm 3kHz |
| Main unit weight: | 1.35kg |
| Adapter weight: | 1.2kg |
| Dimensions: | 237mm \times 206mm \times 95mm |
| Tank capacity: | 400 ml |



Medical equipment for veterinarians

22 rue Charles de Gaulle
22100 Lanvallay - France
Tél. +33 (0)2 96 85 86 76

export@manomedical.com

www.manomedical.com