

Coag&Cut M200

**200W Monopolar & Bipolar
Electrosurgery unit**

Perfect tool for intensive daily use



High Performance





- **5" LCD Touch Screen :**
Multilingual interface available in English, French and Spanish
- **Included accessories :**
 - ▶ A monopolar handswitch handpiece
 - ▶ Five electrodes (straight & angled blades, ball, loop and needle)
 - ▶ A rubber patient return electrode
 - ▶ A bipolar handpiece
 - ▶ Double footswitch integrating a mode switching control button
- **Dimensions & Weight :**
 - ▶ 38x31x12 cm
 - ▶ 6 kg

6SENSE Technology

For more sensibility

The real SPRAY Coagulation

Homogeneous, Superficial & Faster coagulation



Smart surgical assist

Program Manager allows smart access to programs & settings

Self Diagnosis

Quick identification of errors and immediate solution

▶ Maximum Power (Watts) according to the selected mode :

Monopolar Cut	Low	200W @300Ω
	Pure	200W @300Ω
	Blend	120W @300Ω
Bipolar Cut	Macro	50W @100Ω
Monopolar Coag	Soft	80W @500Ω
	Fulgurate	80W @500Ω
	Spray	80W @500Ω
Bipolar Coag	Micro	50W @100Ω
	Standard	50W @100Ω
	Force	50W @100Ω

▶ Input Power :

Operating Voltage: 180-265Volt AC, 50-60 Hz
Max Power input : 1000 VA



22 rue Charles de Gaulle
22100 Lanvallay - France
Tel. +33 (0)2 96 85 86 76

export@manomedical.com

www.manomedical.com