

Blood Pressure

petMAP *graphic II*

10:42

Cat Tail

115 SYS

79 DIA

87 MAP

113 HR

3 Min BPcycle OFF LoAlarm OFF HiAlarm

***** Time to start new BP *****

40	115 / 079	(089)	113
39	157 / 082	(107)	102
39	098 / 074	(081)	106
38	097 / 080	(085)	116
38	171 / 144	(153)	113
08	103 / 083	(089)	116
07	114 / 074	(086)	116
06	189 / 148	(161)	121
06	133 / 060	(079)	120
05	100 / 067	(079)	118

Other Species CuffSite BP Go/Stop

petMAP
graphic II



petMAP *graphic II*

Blood pressure device

Designed for all your BP/HR measurement needs

PetMAP™ graphic II has been specially developed for a veterinary use and offers the PPO technology, which allows to measure:

- Systolic and diastolic blood pressure
- Mean arterial pressure (MAP)
- Heart rate

It comes with 7 cuffs of different sizes.



- ▶ Display of systolic, diastolic, mean arterial pressure (MAP) along with heart rate
- ▶ Correlation to intra-arterial readings
- ▶ Easy to use «no menu» touchscreen display
- ▶ Real time display of oscillometric «envelope»
- ▶ Selectable BP cycle times for measurement reminders
- ▶ Graphic display of last 45 readings and NSV
- ▶ Tabular display of readings
- ▶ Alarm set for systolic blood pressure
- ▶ BP cycle time set for automatic reminders to start another BP measurement
- ▶ User supplied SD memory card records all BP/HR values and «power off» images

Options:

- Large cuffs: 6.5 - 8 - 10 and 13 cm
- Protective cover
- Pole mount
- Table mount
- AC/DC adapter

Technical data:

▶ Systolic BP:	40 - 265 mmHg
▶ MAP:	30 - 200 mmHg
▶ Diastolic:	25 - 195 mmHg
▶ Heart rate:	40 -240 bpm
▶ Temperature:	29°C - 46°C
▶ Touchscreen	4.3" (10.9 cm)
▶ Dimensions and weight:	15 x 8 x 3 cm - 400 gr (piles comprises)
▶ Batteries:	4 AA Alkaline
▶ Included cuffs:	2 - 2.5 - 3 - 3.5 - 4 - 4.5 and 5.5 cm



22 rue Charles de Gaulle
22100 Lanvallay - France
Tel. +33 (0)2 96 85 86 76

export@manomedical.com

www.manomedical.com